
Compact and Light-Weight Automated Semen Analysis Platform Using Lensfree on-Chip Microscopy

Author: Ting-Wei Su,† Anthony Erlinger,† Derek Tseng,† and Aydogan Ozcan*,†,‡


Electrical Engineering Department, University of California, Los Angeles, California, and California NanoSystems
Institute, University of California, Los Angeles, California

Presenter: DUCK-HA HWANG













Birth Rates, 1920-2010

Births per 1,000 women ages 15-44



Source: Statistics calculated using data obtained from the National Center for Health Statistics and Heuser (1976), available here

PEW RESEARCH CENTER





UCLA

HOW?

Analysis Sperm's fertility

 Digital subtraction of these consecutive lensfree frames, followed by appropriate processing of the reconstructed images

Density


- Count sperm on the image. (High density, High fertility)

Motility


- Count moving sperms, Count stationary sperms. (High motility, High fertility)

Speed

 Calculate speed using the dynamic trajectories of motile sperms(High speed, High fertility)




Recorded Hologram **Zoom In UCLA**




Reconstruct Image (Computational Process)

Phase

Amplitude





UCLA




Sperms' Head


Phase Image





Sperms' 2D Trajectory










Chamber



Cell VU 600 Depth 20 um

Size of Human Sperm's Head head 5.1 μ m by 3.1 μ m and a tail 50 μ m long

HOW?

Subtract Background

(ex: First Hologram Image – Second Hologram Image)



Reconstruct Image



Find moving objects




Compare each pixel




Record trajectory






Analysis Sperms



Count (number of sperm): 233 Median (speed): 16.4004 um/sec Mean (speed): 21.9802 um/sec



Motility

Stationary sperms

Average hologram image


(Combine each frame and divided by number)



leave stationary ones





Count stationary sperms




Moving sperms

Stationary sperms





Analysis Fertility of sperms





Commercialization

- More widespread Lens-free Microscopes
 - -Cost effective
 - -Easy user interface
 - -Fast analysis
 - -Simple product
 - -Design

Conclusion

END

